



ProGuide
implant excellence, guided with care

ProGuide Process



1- You will need to send us a panoramic, a plaster model or a digital impression* and a submission form.

** Digital Impression is sent through Trios direct connection to "ProGuide" lab connection by adding the account under proguide.prodent@gmail.com or STL files are sent by email/we transfer to proguide.prodent@gmail.com.*

2- ProGuide team will provide the adequate process to follow for the CBCT in 24hrs.

3- A surgical report and a quotation will be sent to your approval after 72 hours of receiving the plaster model/digital impression and CBCT.

4- In case any modification is required, the revised surgical report will be sent within 24 hours.

5- Approval is taken into account when you reply to the email or sign hard copies of both aforementioned documents.

6- Upon approval, the Surgical Guide will be manufactured within 72 hours.

7- The drilling protocol and the surgical report will be sent with the guide.

8- Our technical team or customer service will contact you before the Surgery date.

ProGuide Submission Form

| PATIENT NAME | JAW | IMPLANT BRAND | # OF IMPLANTS | SITE NUMBER | SURGERY DATE |
|--------------|-----|---------------|---------------|-------------|--------------|
| | | | | | |
| | | | | | |
| | | | | | |

Guided Surgical Kit: Yes No

Comments: